# IVDK WinAlldatNet

# Administrator's Guide

Instructions for converting the database system to Microsoft SQL Server

5th revised edition 2024-05-22 Autors: CM, RA

# 0.0 Preface

The conversion to another database system should only be done in consultation with and together with an IVDK IT employee.

# The abbreviation CWS = clinical workplace systems is used several times

# **1.0 Requirements**

- Installed SQL Server (with costs, NOT provided by IVDK!) or
- SQL Server Express (free of charge). The latter does not require a server, but can also be installed on a normal Windows PC!
- A WinAlldatNet installation with access to the Access database used so far.
- Windows user profile with sufficient permissions to create a database on the SQL Server intended for WinAlldatNet.
- Access to an installed SQL Server Management Studio (SSMS) to set user permissions on the SQL Server.

# **1.1 Terminology**

Instance name (assigned during SQL Server installation)
 When using the default instance, an SQL Server is addressed by the name of the Windows Server, alternatively the IP address of the Windows Server can be used. If several SQL Servers are running on one server, they are distinguished by the instance name. For the database connection, the instance name is appended to the server name with \: Server name or IP address of the server instance name.

For SQL Server Express, a non-changeable instance name is usually assigned. Instance names other than the default instance reduce portability, especially if the Integration or Analysis services of SQL Server (Not included in Express and Standard!) are used!

• IS (Integrated Security)

With Integrated Security, the current Windows user login is used as the login to the SQL Server database. The specification of a user name, as well as the associated password are omitted here. The SQL server still has the older procedure, with which the user names together with password are deposited on the SQL server (mixed logon, only with the installation adjustable!), this older procedure is however more insecure, since with the data base logon the passwords are transmitted unencrypted!

The WinAlldatNet program uses the Integrated Security by default to logon to the SQL Server.

• Database

You can create several databases on one SQL Server. These databases are managed by the master database, are physically separated (My Documents) and also have separate access permissions. For administrators who know only Oracle databases: Oracle does not know databases. There the visibility and the access permissions of the objects are defined by schemas alone. Access across databases (if authorization is given!) is done via database name.schema.object (object = table, view etc.). The specification of the schema is omitted if

schema dbo.

• Schema

SQL Server uses dbo (DataBase Owner) as default schema to access the tables and other database objects. Users or groups can be assigned a schema other than dbo. Different schemas use a separate table for each table name, which is separated from the others and has the access permissions of the respective schema. For example, the schemes dbo, department1, department2 and department3 would result in the separate tables dbo.invoices, department1.invoices, department2.invoices and department3.invoices for the table invoices. In case of schema foreign access, the schema must be preceded by a . separated from the object (object = table, view, etc.) and requires the appropriate access permissions. WinAlldatNet uses only the default schema dbo. Another schema is not provided and is not supported by the program!

- Alias (Assignment via SQL Server(Configuration)Manager) Instead of IP address of the server instance name, the access to a SQL server can be simplified by assigning an alias name e.g. winalldatnet. If frequent restructuring of the servers takes place, it makes sense to assign an alias name, since the access name (alias) does not change when the database is moved to another server and thus does not need to be changed for the accessing CWS. Unfortunately, every now and then the network administrators manage to block the name resolution of the alias names of the SQL Server via the network. Therefore, when using an alias name, always check the alias name resolution first in case of connection problems!
- Localhost localhost are assigned alias names for access to locally located SQL Servers. The name abbreviation applies only to the server name or IP address of the server. Instance names must still be appended with \InstanceName.
   Can only be used if the SQL Server and the WinAlldatNet program are located on the same server/PC (SQL Server Express).

# 1.2 Note for SQL Server installation

SQL Server uses several services that require corresponding administrative user accounts. If possible, use local accounts of the server / CWS. If a domain-only account is used, and it no longer exists, the SQL Server can no longer be administered and must be reinstalled!

# 2.0 Preparations

# 2.1 Global folder

A global folder on a network drive is required under which the WinAlldatNet program can store shared files. All WinAlldatNet users must have read and write access to this folder and must be able to create new file folders. This global folder serves among other things to contain the winalldatnetpwd.ini for the user login, the log files for program errors and the backups of the DKG updates. For simplicity, it is recommended to use the existing Access database folder for this purpose, which previously fulfilled this task.

# 2.2 Domain group for WinAlldatNet users (optional)

A domain group in which the WinAlldatNet 2 users are entered has proven to be useful. If a user changes, no changes have to be made to the SQL Server.

# 2.3. Configuration of SQL Server

Open the SQL Server(Configuration)Manager on the server/PC with the SQL Server. This is no longer linked with some SQL Server versions, but must be called manually under C:\Windows\System32\SQLServerManagerXX.msc. The XX in the file name stands for the SQL Server version number. If more than one version is available, the newest version must be started.

Datei Aktion Ansicht ?	
SQL Server-Konfigurations-Manager (Lokal)       Name       Status       Status       Anmelden als       Proze         SQL Server-Dienste       SQL Server-Netzwerkkonfiguration (32 Bit SQL Server Analysis Services (SQLSERVER12)       Wird ausgeführt       Automatisch       NT AUTHORITY/NE       11076         SQL Server-Netzwerkkonfiguration       SQL Server Analysis Services (SQLSERVER12)       Wird ausgeführt       Automatisch       NT Service/MSOLA       5000         SQL Server (SQLSERVER12)       Wird ausgeführt       Automatisch       NT Service/MSOLA       5000         SQL Server (SQLSERVER12)       Wird ausgeführt       Manuell       NT AUTHORITY/LO       3440         SQL Server (MSSQLSERVER)       Wird ausgeführt       Manuell       NT Service/MSSQL       5304         SQL Server Reporting Services (MSSQLSERVER)       Wird ausgeführt       Automatisch       LocalSystem       3908         SQL Server Reporting Services (SQLSERVER)       Wird ausgeführt       Automatisch       LocalSystem       0         SQL Server Reporting Services (SQLSERVER)       Beendet       Andere (Boot, Syste       LocalSystem       0         SQL Server Reporting Servi	ss-ID
SQL Server-Agent (MSSQLSERVER) Wird ausgeführt Automatisch LocalSystem 3824	
	>

Under "SQL Server Services", in addition to the actual SQL Server, the services "SQL Server Browser" and the "SQL Server Agent" associated with the SQL Server should run.

🧱 Sql Server Configuration Manager				-	×
Datei         Aktion         Ansicht         ?					
<ul> <li>SQL Server-Konfigurations-Manager (Lokal)</li> <li>SQL Server-Dienste</li> <li>SQL Server-Netzwerkkonfiguration (32 Bir</li> <li>SQL Native Client 11.0-Konfiguration (32-</li> <li>SQL Server-Netzwerkkonfiguration</li> <li>Protokolle für 'MSSQLSERVER'</li> <li>Protokolle für 'SQLSERVER2'</li> <li>SQL Native Client 11.0-Konfiguration</li> </ul>	Protokollname Shared Memory Named Pipes TCP/IP	Status Aktiviert Deaktiviert Aktiviert			

In the "SQL Server network configuration" or "SQL Server network configuration (32-bit)", check whether the TCP/IP protocol is enabled for the SQL Server in question. Alternatively, the obsolete "Named Pipes" protocol can be deactivated. A database connection via Named Pipes is about 50% slower than via TCP/IP.

For the TCP/IP protocol, the properties must be called up (right mouse button on the entry). In the last entry (IPAll), under "Dynamic TCP ports", the entry is to be emptied and in the "TCP port" field below it, 1433 is to be entered. This procedure facilitates the setting of the firewall parameters for WinAlldatNet installations.

	IP-Adresse	fe80::5efe:10.2.3.5%5	^
	TCP-Port	1433	
Ξ	IP6		
	Aktiv	Ja	
	Aktiviert	Nein	
	Dynamische TCP-Ports		
	IP-Adresse	2001:0:5ef5:79fd:2496:c56:a7b9:f;	
	TCP-Port	1433	
-	IP7		
	Aktiv	Ja	
	Aktiviert	Nein	
	Dynamische TCP-Ports		
	IP-Adresse	fe80::2496:c56:a7b9:fa57%6	
_	TCP-Port	1433	
	IPAII		
	Dynamische TCP-Ports		
	TCP-Port	1433	
<b>Al</b> Gi	<b>ctiv</b> bt an, ob die ausgewählte IP	-Adresse aktiv ist.	

Proceed similarly with the "Native-Client Configuration" or "Native-Client Configuration (32 Bit)".



Also there you have to check under properties of the TCP/IP protocol if the SQL Server standard port 1433 is entered.

Eigenschaften von TCP/IP		?	×
Protokoll			
□ Allgemein			
Aktiviert	Ja		-
Erhalten	30000		
Erhaltungsintervall	1000		
Standardport	1433		
<b>Aktiviert</b> TCP/IP-Protokoll für diesen (	Client aktivieren oder (	leaktivieren	
ОК	Abbrechen Über	rnehmen	Hilfe

# 2.4 Configuration firewall for the SQL Server

On the server/PC on which the SQL Server with the WinAlldatNet databases is to run, the following ports must be released in the firewall: TCP port 1433 (SQL Server connection) TCP/UDP port 1434 (SQL Server browser)

# 2.5 Configuration firewall for SQL Server Express on normal Windows PC

As of Windows Vista, the Windows operating systems have their own operating system firewall. TCP port 1433 must be enabled there.

Zugelassene Apps			- 0	
→ 🔹 🛧 🝻 « Windows Defender Fir > Zugelassene Apps	∿ Ç			
Communikation von Apps durch die Windows Defender licken Sie zum Hinzufügen, Ändern oder Entfernen zugelassener Apps ur Velche Risiken bestehen, wenn die Kommunikation einer App zugelassen vird?	Firewall : nd Ports auf	zulas "Einst Eins'	sen ellungen änd tellungen än	lern' derr
Zugelassene Apps und Features:				
Name	P	rivat	Öffentlich	^
Skyrim Herbarium		✓		
SNMP-Trap				
SonicWALL Mobile Connect		✓	$\checkmark$	
SonicWALL.MobileConnect		✓	✓	
✓ Spielleiste		✓	$\checkmark$	
Sprachausgabe-Schnellstart		✓	◄	
SQLServer				
Ssh Server		✓		
☑ Sway		✓	✓	
✓ System Test		✓		
✓ Take a Test		✓	$\checkmark$	
Teamviewer Remote Control Application			✓	~
	Details	s	Entferne	n
	A	ndere	App zulasse	n
	OK		Abbrech	ien

Port bearbeiten X							
Verwenden Sie diese Einstellungen, um einen Port durch die Windows Defender Firewall zu öffnen. Wählen Sie einen sinnvollen Namen aus, den Sie sich merken können. Informationen zu der Portnummer und zu dem Protokoll finden Sie in der Dokumentation der App bzw. des Diensts, die bzw. den Sie verwenden möchten.							
Name:	þQLServer						
Portnummer:	1433						
Protokoll:	● TCP						
Welche Risiken bestehen beim Öffnen eines Ports?							
Sie können auswählen, für welche Netzwerktypen der Port geöffnet ist.							
Netzwerktypen OK Abbrechen							

Alternatively, the following Windows shell script can be executed for the setting:

@netsh advfirewall firewall add rule name = SQLPort dir = in protocol = tcp \action = allow localport = 1433 remoteip = localsubnet profile = DOMAIN

# 2.6. Configuration Firewall for the WinAlldatNet Clients

The CWS with WinAlldatNet require an unlocked TCP port 1433.

#### 3.0 Immediate preparations before the changeover

- All WinAlldatNet users have to stop data input and close all open WinAlldatNet programs.
- On the PC with the WinAlldatNet installation from which the conversion takes place, the following line must be added to the winalldatnet.ini configuration file:
   *-localimportexportadmin* Note: The entry must be in a separate line!

#### 4.0 Performing the database conversion

Log on to the PC with WinAlldatNet installation with the Windows user profile, which has sufficient permissions to create a database on the SQL server intended for WinAlldatNet. Now start WinAlldatNet with the user Chef (password IVDK if it has not been changed). Now open the submenu "Export data" in the menu "Utilities" and select the database system into which the currently used database should be converted.

WinAlldat									
Daily routine	Adı	ministration	Languages	Documents	Extras	Help	)		
		Test series n	Test series management						
Nev		Deactivate t	est series		Exam	ination date	03.05.2021	<ul> <li>Station</li> </ul>	
Patient	82	Data backup	p	•					Advanced Search
E	$\mathbf{X}$	Clean and r	epair database	1				~	
N	$\times$	Delete curre	ent case / patie	ent	anofor d	lata to	the IVDK		Filter
		Test substar	nce list test seri	ies					Theor
Ana		Test substar	nce list with filt	ter		Tes	st series in test	plan	O Type 1
		Test series li	st						U 17F-1
Add pat		List of multi	iple-tested sub	ostances	-				Test
Add pa		Relationship	o testmethod-	relevance					method
Select DK0		Additional o	questions						
	۵.	Access auth	orization						
Text b		Import data	1						
		Export data		•	-:	> SQLit	te		
Readings		Extract files	(ZIP, 7-Zip)		-1	> SQL S	Server		
	_				-:	> Acce	SS		
			Name		-:	> text			
					-)	> text (	tables separa	ted)	
Pelaurana dialag				-)	> Text (	(HFG only)			
Relevance dialog					-)	> Text (	(empty HFG a	and HFG tables	s only)
Final doo		ntation			D	KG Upo	date Access -	> SQLite	
Printed output	ıt	File export	1	Fest protocol form		Releva	ance form	Test p	rotocol

Now enter the IP address or the server name of the SQL server. If additional instance names were assigned to the SQL server, the instance name must be appended to the IP address or the server name as \instance name.

🔒 Server log	in —		×
Server:	IP or server name	~	0
Password:			
	Integrated security		_
ОК		Car	ncel .:

Alternatively, you can use the information icon to search for the server in the network. This search may take some time. The search process is finished as soon as the combo box can be opened or opens automatically with newer WinAlldatNet versions.

The checkbox Integrated Security uses the current Windows user login as login to the SQL Server database. The specification of a user name and the corresponding password is not necessary.

🧧 Copy database data	1			-	×
Table	Read	Write	Status		
		C	opy data		

To start the conversion process press the "Copy data" button. Non-existing databases will be created, if they do not exist. Before the data transfer, the database update is performed, thus bringing the target database up to date.

During the run, the conversion status is displayed for each table. Successfully converted tables get

the status "Ok". Tables that do not exist in the target database get the status "Table not exists". These tables are not used by WinAlldatNet 2 or are located in an external database. In case of a faulty conversion, the faulty SQL statements are displayed as status text. A double click with the left mouse button shows the error text in a separate text display.

In the case of the TRANSACTION\_LOG table, error messages may be generated by SQL Script in the SQL. This is only the case if logging according to the Swiss Human Research Act (HRA) is switched on in the options.

🥫 Copy database da	ita			- 0	×
Table	Read	Write	Status		^
tblPatblock	7091	7091	ОК		
tblPatblocksubst	25601	25601	ОК		
tblPatchTest	51602	51602	ОК		
tblPatient	13447	13447	ОК		
tblPatientNew	0	0	Table not exists		
tblPatientOld	189	189	ОК		
tblPatient Tmp	0	0	ОК		
tblReaction	32164	32164	ОК		
tblReactionCat	88	88	ОК		
tblReading	102540	102540	ОК		
tblReadingNote	12604	12604	ОК		
tblRelevanceCat	88	88	ОК		
tblSelectionCat	1656	1656	ОК		
tblStationCat	26	26	ОК		
tblSubstance	3478	3478	ОК		~
			Consudate		
			Copy data		
					.:

# 5.0 Final work

The former Access database should now be archived and removed or renamed so that no WinAlldatNet installation can enter new data there.

# 5.1 Setting the database authorizations in SQL server

Now start the SQL Server Management Studio. This can be located on the server/PC with the SQL Server or as a separate installation on another PC. Connect to the SQL Server that contains the WinAlldatNet databases and select the Logins section under Security in the Object Explorer. Note that this is the security of the server and not the security area of the respective databases!



# 5.1.1 Adding a new user / group

If the users or the user group for WinAlldatNet are not yet entered under SQL Server logins, open the selection menu with the right mouse button of the logins and select "New Login...".

🗄 Login - New				-	
Select a page	T Script 🔻 🖓 Help				
<ul> <li>General</li> <li>Server Roles</li> <li>User Mapping</li> <li>Securables</li> <li>Status</li> </ul>	Login name: <u>Windows authentication</u> <u>SQL Server authentication</u> <u>Password</u> : <u>Confirm password</u> :				Search
	Old password: Enforce password policy Enforce password expira User must change passw Mapped to cettificate Mapped to asymmetric key	tion vord at next login		~	
Caracter	Map to Credential			~	Add
Server: CHR11	Mapped Credentials	Credential	Provider		
Connection: CHR11\Meyer_2					
v₩ <u>View connection properties</u>					
Progress					Remo <u>v</u> e
Ready	Default <u>d</u> atabase: Default lan <u>q</u> uage:	master <default></default>		~ ~	
			[	ОК	Cancel

In the newly opened "Login - New" window, click the "Search ..." button to search for the user / user group in the domain path.

Benutzer oder Gruppe auswählen	×
Objekttyp: Benutzer oder Integriertes Sicherheitsprinzipal	Objekttypen
Suchpfad:	
CHR12	Pfade
<u>G</u> eben Sie die zu verwendenden Objektnamen ein ( <u>Beispiele</u> ):	Namen überprüfen
Erweitert OK	Abbrechen

If you want to add a user group, click the "Object types ..." button. Select "Groups" as the object type. This is usually not selected!

Objekttypen	×
Wählen Sie den Objekttyp, den Sie suchen möchten. Objekttypen:	
Andere Objekte       Integrierte Sicherheitsprinzipale       Gruppen       Benutzer	
OK Abbrechen	

Back in the "Select user or group" window, you also have the option under "Paths..." to set the starting point of the search in the domain structure. You can now enter the user name or the name of the user group under "Enter the object names to be used" and click on the "Check names" button. Alternatively, you can open the advanced search ("Advanced..." button).

Benutzer oder G	uppe auswählen		×
Objekttyp:			
Benutzer, Gruppe	oder Integriertes Sicherheitspri	nzipal	Objekttypen
Suchpfad:			Díada
Allgemeine Abfra	gen		Flade
Name:	Beginnt mit V		Spalten
Beschreibung:	Beginnt mit 🗸		Jetzt suchen
Deaktivierte	Konten ende Kennwörter		Beenden
Tage seit der le	tzten Anmeldung: 🛛 🗸		<del>?</del>
Suchergebnisse:			OK Abbrechen
Name	Ordner		

The "Search now..." button displays the available users and groups, depending on the selected object types.

Benutzer oder Gruppe auswä	hlen			×
Objekttyp:				
Benutzer, Gruppe oder Integrie	rtes Sicherheitsprin:	zipal	(	Objekttypen
Suchpfad:				
CHR12				Pfade
Allgemeine Abfragen				
Name: Beginnt mit	$\sim$			Spalten
Beschreibung: Beginnt mit	$\sim$			Jetzt suchen
Deaktivierte Konten				Beenden
Nicht-ablaufende Kennw	örter			
	_			<b>67</b>
I age seit der letzten Anmeld	ung: 🗸 🗸			7 p
Suchergebnisse:			OK	Abbrechen
Name	Ordner			^
Administrator	CHR12			
Administratoren	CHR12			
ANONYMOUS-ANMELDUN	G			
BATCH				
Benutzer	CHR12			
Scmeyer	CHR12			
DefaultAccount	CHR12			
M DIALUP				~

If you have one or more domains, you must open the corresponding branches to get to the desired user or group. Depending on the permissions, a user query may appear to open such a branch. When you have selected the appropriate user or group, click "Ok".

Benutzer oder Gruppe auswählen	×
Objekttyp: Benutzer oder Integriertes Sicherheitsprinzipal	Objekttypen
Suchpfad:	
CHR11	Pfade
<u>G</u> eben Sie die zu verwendenden Objektnamen ein ( <u>Beispiele</u> ):	
X220-CM\cmeyer	Namen überprüfen
Erweitert OK	Abbrechen

The selected user or user group is now entered and the creation can be completed with "OK" and you are back in the "Login - New" window to set the permissions.

#### 5.1.2 Set the user permissions

If you have skipped 5.1.1 because the user or user group is already entered under SQL Server Logins, open the "Logins" branch selected under 5.1 in the Object Explorer and select the desired user or group by double-clicking on it.

Now switch from to the "User assignment".

E Login Properties - CHR11\	Meyer			- C	I X
Select a page	C Script	🔻 🙆 Help			
🖉 General	a) conpr	•			
Server Roles		and the state to state			
Securables	Users ma	ppe <u>a</u> to this login:	11	D.C. h.C.I	
✗ Status	мар	Database	User	Derault Schema	
-		master			
		model			
		msdb			
		ReportServer			
		ReportServerTempDB			
		tempdb			
		limebase			
		WinAlldatIvdk	CHR11\Meyer		
	$\sim$	Win Aldat Ivdk Report	CHR11\Meyer		
Connection					
	Guest	account enabled for: WinA	ldatlvdk		
CHR11	D.1.1				
Connections	Database	e role membership for: WinAl	dativdk		
CHR11\Meyer 2	db_a	ccessadmin			
		ackupoperator atareader			
TT <u>wew connection properties</u>	dbd	atawriter			
	db_d	dladmin			
		enydatareader enydatawriter			
Progress	∏ db_0	wher			
A Ready	db_s	ecurityadmin			
neauy	🔽 public	2			
, d b .					
				OK	Cancel

There the databases "WinAlldatIvdk" and "WinAlldatIvdkReport" are to be assigned. The WinAlldatNet program needs the dbo schema as default schema, because database updates are imported via the WinAlldatNet program.

In order to add the default schema, press the respective button "..." belonging to the database.

最 Select Schema	×
Select these object types:	
Schemas	Object Types
Enter the object names to select ( <u>examples</u> ):	
	Check Names
	Browse
ОК С	Cancel Help

"Schemas" is already selected as the object type, so you can enter the schema dbo under "Enter the object names to be used". You must then verify the name with "Verify name". Alternatively, you can select the schema via "Browse...".

🖧 Browse	e for Objects		×
13 objects	were found matching the types you selected.		
Matching (	bjects:	Turne	
<b>—</b> 4	[db_depydatawriter]	Sche	
	[db_owner]	Sche	-
	[db_securityadmin]	Sche	_
	[dbo]	Sche	
	[auest]	Sche	
	[INFORMATION SCHEMA]	Sche	
<u> </u>	lough	Sobo	4
	OK Cancel	Help	

Finally, for both databases, under "Membership in database role for: Database name", select "db\_owner".

. Login Properties - CHR11\M	leyer			_	
Select a page		🔻 🙆 Help			
🔑 General	<ul> <li>South</li> </ul>				
Server Roles					
Securables	Users ma	pped to this login:		<b>D</b> ( ) <b>D</b> (	
Status	Map	Database	User	Default Schema	
		master			
		model			
		msdb			
		ReportServer			
		ReportServerTempDB			
		tempdb			
		Timebase			
		WinAlldatIvdk	CHR11\Meyer	dbo	
	$\checkmark$	Win Alldat lvdk Report	CHR11\Meyer	dbo	
Connection					
Connection	Guest	account enabled for: WinAl	datlvdk		
Server: CHR11	Database	role membership for: WinAll	datlvdk		
Connection:	db_a	ccessadmin			
CHR11\Meyer_2	db_bi	ackupoperator			
<u> </u>	db_da	atareader			
		tawnter Iladmin			
	dbd	enydatareader			
_	db_de	enydatawriter			
Progress	db_o	wner			
Ready	db_se ✓ public	ecuntyadmin :			
				ОК	Cancel

All other required settings correspond to the default settings and do not need to be changed. With "OK" the changes are completed.

Note: It is possible to use the local SQLite database winalldatnetrep.db3 instead of the WinAlldatIvdkReport database for "Reports and Lists". The settings for the WinAlldatIvdkReport database would be omitted in this case!

#### **5.1.3** Alternative settings of user permissions

Due to regulations of individual data protection officers, it may not be possible to set the schema "dbo" and the database role "db\_owner". In this case the schema remains empty. Two variants are possible as database role entries:

Variant 1: Read/write authorization and authorization to change tables

For this, instead of the database role "db\_owner", the database roles "db\_Datareader", "db Datawriter" and "db dlladmin" must be selected.

🗄 Login - New				—	
Select a page		🗙 🙆 Help			
🖉 General	2 John	. Top			
<ul> <li>Server Roles</li> <li>User Mapping</li> </ul>	Users ma	oped to this login:			
<ul> <li>Securables</li> <li>Status</li> </ul>	Map	Database	User	Default Schema	
		model			
		msdb			
		Report Server			
		Report Server Temp DB			
		tempdb			
		Timebase			
		WinAlldatlvdk			
		Win Alldat lvdk Report			
Connection					
Connection	Guest	account enabled for: WinAll	datlvdk		
CHR11	Database	role membership for: WinAll	latlvdk		
Connection:	db_a	ccessadmin			
CHRITWeyer_2	db_bi	ackupoperator			
View connection properties	⊡ db_da	atareader Hawriter			
		lladmin			
	db_de	enydatareader			
Progress		enydatawriter vper			
Dead.	db_se	curityadmin			
neauy	✓ public	:			
				ОК	Cancel

In this variant there are hardly any restrictions for the operation of the program. Only database backups cannot be created from within the program.

Variant 2: Only read/write authorization without authorization to change tables.

For this, instead of the database role "db\_owner", the database roles "db\_datareader", and "db\_datawriter" must be selected.

🗄 Login - New			- 0	×
Select a page	🖵 Script 🔻 <table-cell> Help</table-cell>			
<ul> <li>General</li> <li>Server Roles</li> <li>User Mapping</li> </ul>	Users mapped to this login	:		
<ul> <li>Securables</li> <li>Status</li> </ul>	Map Database master model msdb ReportServer ReportServer tempdb Timebase WinAlldatIvdk	iempDB Report	Default Schema	
Connection Server:	Guest account enabled	for: WinAlldatlvdk		
CHR11 Connection: CHR11\Meyer_2 V Vew connection properties	Database role membership db_accessadmin db_backupoperator db_datareader db_datareader db_ddladmin db_denydatareader db_denydatareader db_denydatareader	for: WinAlldatlvdk		
Ready	☐ db_securityadmin ☐ db_securityadmin ☑ public			
			ОК С	Cancel

For this variant the WinAlldatNet program must be additionally called with the parameter -paranoia. This parameter locks program areas which require authorizations to change database tables and issues a corresponding message. These program areas must then be executed by a user with appropriate administrative database authorizations. Affected program functions are the DKG test series update, the database update after importing a new program version and the creation of database backups from within the program.

This parameter prevents users from executing program functions that make changes to database tables. In this case, a corresponding message is issued. The DKG test series update, database update after program update and database backup cannot be executed by the normal users from the program. For this purpose, an administrative user with corresponding authorizations on the database is absolutely necessary!

# 5.2 Setting up the clients Accessing the databases via winalldatnet.ini

#### 5.2.1 Introduction to the structure of the winalldatnet.ini configuration file

WinAlldatNet no longer uses the ODBC interface of its predecessors WinAlldat IVDK and WinAlldat ESSCA for the database connection. For compatibility reasons, however, the file path to the database is still read from the ODBC interface first and used as connection to the Access database. Under WinAlldatNet 2 the database parameters of WinAlldatNet.exe are passed as command line parameters. The available parameters are displayed with WinAlldatNet.exe -? To simplify maintenance for multi-user installations, command line parameters can also be passed using the winalldatnet.ini text file. It should be noted that all parameters used and their possible sub-parameters must be in separate lines.

Example:

winAlldatNet.exe -log -dbpath d: \ winalldat \ winalldativdk.mdb -jet is entered in winalldatnet.ini as follows:

-log
-dbpath
d:\winalldat\winalldativdk.mdb
-jet

If the same parameter occurs both as a command line parameter and in winalldatnet.ini, the command line parameter is overwritten by the winalldatnet.ini entry. The winalldatnet.ini is not case sensitive. All - parameter beginnings can also be replaced by / (Ex: -log can also be written as /log).

#### 5.2.2 Client access via Integrated Security (SQL server)

The user/domain user logged on under Windows is used here to log on to the SQL Server.

Content of winalldatnet.ini:

-sqlserver -server IP of the server or server name\instance name (\instance name is omitted if not available) -globalpath Drive letter:\file path global folder\

Alternatively for drive letter: \file path global folder\ you can also use \\servername\file path global folder\

The global folder on a network drive is needed so that the WinAlldatNet program can store shared files. All WinAlldatNet users must have read and write permissions there and must be able to create new file folders. This global folder is used among other things to hold the winalldatnetpwd.ini which contains the hash values for the WinAlldatNet user logins. For simplification the already existing network folder on which the Access database was located can be used.

# 5.2.3 Client access with restricted user permissions

# 5.2.3.1 Operating system

Client access for installation environments without write file permissions to the WinAlldatNet program folder (also applies to installations under c:\Programs!)

#### -readonlyinstall

This parameter causes all files that are changed by the program during operation to be swapped out to the global folder. The parameter

-globalpath with drive path as sub-parameter is mandatory additionally!

# 5.2.4 Other useful winalldatnet.ini parameters

-log

WinAlldatNet logs all exception errors. With the parameter -log additionally all database accesses with errors (read and write) are also logged.

#### -update

WinAlldatNet checks after the program call whether a new version is available. If a new maintenance update is available, it will be downloaded via HTTP protocol from the IVDK server and installed.

If the HTTP protocol on the clinical workstation systems (CWS) is blocked by a firewall, the intranet update can optionally be activated by specifying the path (windowshare) of an intranet server.

#### -update

Variant 1 (classic):	Drive letter:\file path intranet folder
Variant 2 (UNC):	\\ servername\file path intranet folder

The following screenshot shows the classic setting (variant 1) as an example:

🔀 Options	×	ζ
Other options Design Printouts Keyb	oard, anamnesis and text modules Data transfer KIS	
Data transfer		
○ SMTP (e-mail) ○ FTP	Other transfer mode	
Clinic name	×	
Transfer folder	C:\winalldatnet\DataExport	
14 👻	Number of days between data transfers (IVDK recommendation: 14 days)	
SMTP		
Name of the SMTP server		
E-mail address sender		
Name of the sender		
Port (25, 465, 587)	25 文 Timeout / s 30000	
Login name for SMTP server		
SSL / TLS encryption	Copy to sender	
Send emails to folder		
Internet update		
Proxy Name / IP	Port 80	
FTP		
FTP clinic folder		
Proxy Name / IP	Port 80 Encrypted transmission (SFTP)	
Intranet Update		
Global Update Folder	D:\winalldatnet\Update	
Save	Cancel	

The current WinalldatNet maintenance update must be copied manually into the intranet folder. The WinAlldatNet update consists of the 2 files *version.txt* and *winalldatnet.zip* and can be downloaded from the IVDK server at *https://winalldatnet.ivdk.org/update/*.

*-x64* 

This parameter is mandatory if a 64-bit Microsoft Office is installed!

# -localreport

When using SQL Server as database, a separate SQL Server database "WinAlldatIvdkReport" is used for reports and lists. With this parameter the local winalldatnetrep.db3 is still used for reports and lists.

# 5.2.5 winalldatnet.ini Parameters for Administrators

#### -dbtest

This parameter starts a test window with checking of possible database connections. Use this parameter on CWS where WinAlldatNet cannot be started because of not found or faulty database. Access errors can be limited in this way. From the test window you can also call the system information built into WinAlldatNet.

#### -localimportexportadmin

This parameter requires the login with the user "Chef"! Unlocks the submenus "Import data" and "Export data" in the "Utilities" menu.

#### -localdbadmin

This parameter requires the login with the user "Chef"!

Unlocks the submenus "DB Editor" and "DB Test" in the "Utilities" menu. Under "DB Test" is the same test window for checking the database connections as when the program is called with the -dbtest parameter. The "DB Editor" is an integrated SQL editor for editing the databases used by WinAlldatNet via SQL commands.

#### 6.0 Appendix

# 6.1 Differences in the database structure Access to SQL Server

When using SQLite as database, the Access file Winalldativdk.mdb, or for ESSCA version Winalldatessca.mdb are replaced by the SQL Server database WinAlldatIvdk. In this case, the tblCaseEssca table included in the IVDK version was replaced by the extended tblCaseEssca table from the ESSCA version.

The SQLite database winalldatnetrep.db3 used for reports and lists has been dropped \*). The tables contained there are located in the SQL Server database WinAlldatIvdkReport. The SQLite ODBC entry needed to access the winalldativdk.db3 database is replaced by a corresponding ODBC entry to the WinAlldatIvdkReport SQL Server database.

\*) Using the -localreport parameter, the locally located winalldatnetrep.db3 can still be used.